

β -Galactosidase Activity Assay

You will need to determine the optimal cell culture conditions for your β -Gal assay. Typically, overnight cultures are diluted 1:50 – 1:100, and then grown to mid-log or early stationary phase. (Overnight cultures can also be used directly for β -Gal assays in some cases.)

Z buffer

To make 500 mL:	Final concentrations:
150 mL 0.2 M Na ₂ HPO ₄	0.06 M Na ₂ HPO ₄
100 mL 0.2 M NaH ₂ PO ₄	0.04 M NaH ₂ PO ₄
5 mL 1 M KCl	0.01 M KCl
5 mL 0.1 M MgSO ₄	0.001 M MgSO ₄

pH should be 7.0; filter sterilize

Just before using, add β -mercaptoethanol (IN THE HOOD) (35 μ L per 10 mL Z buffer).

1. Determine the A₆₀₀ of the cultures
2. In a 2 mL tube, **mix** 0.1 mL cell culture with 0.9 mL Z buffer, 100 μ L chloroform and 50 μ L 0.1% SDS. **Vortex** the solution for 10 sec, then let the tubes **rest** at RT.
3. **Add** 0.2 mL of *o*-nitrophenyl- β -D-galactopyranoside (ONPG) in Z buffer (4 mg/mL) to the lysed culture prepared in step 2. When a yellow color develops, **stop** the reaction by adding 0.3 mL 1 M Na₂CO₃. **Record** the time at which the yellow color stops.
4. Briefly **spin** down the tubes (~1 min full speed; or even 10 seconds in picofuge)
5. Being careful to only take the aqueous (top) phase*, **remove** 1 mL of the assay solution from the tube and **transfer** to a cuvette. **Read** the A₅₅₀ and A₄₂₀ of the reaction samples. Units of β -galactosidase activity =

$$1000 \times \frac{A_{420} - (1.75 \times A_{550})}{t \times 0.1 \times A_{600}}$$

where t = time in minutes.

*The chloroform will mess up the plastic cuvettes

Optional:

Use PopCulture to lyse cells (EMD Biosciences)

Do in 96 well with 140 μ L Z buffer, 10 μ L cell lysate and 10 μ L ONPG. Measure A₄₂₀ for 20-40 minutes at 30-second intervals.