

Phase Lock Gel RNA Preparation

General Notes

- Always practice RNase-Free technique when working with RNA
- Use RNA-only pipet tips and tubes
- Wipe off counter tops and pipetmen with RNaseZap
- Filter-sterilize solutions prior to use
- Use either DEPC-H₂O or MQ-H₂O
- Keep RNA cold (4 °C or lower)

Lysis solution (8x) [make ~ 10 mL at a time]

- 320 mM NaAcetate (pH 5.0)
- 8% SDS
- 16 mM EDTA

Protocol:

- Pellet** bacteria cells by spinning at 4 °C, 5,000 -6,000 x g, for 5 minutes.
 - If growing bacteria in M9 media, pellet *at least* 2 mLs.
 - If growing bacteria in LB media, you can skip this step and just take 500 µL of culture directly to Step 4.
- Remove** supernatant, being careful not to lose pellet
- Resuspend** pellet in 500 µL 1x M9
- Add** 71 µL 8x Lysis Buffer, quickly pipetting up and down to mix. **Immediately add** 570 µL Acid Phenol/Chloroform (Ambion), **vortex**.
- Incubate** 5 minutes at 65 °C with regular **vortexing** (every 50 sec)
- Spin** > 12,000 – 16,000 x g, 10 minutes
- Transfer** the supernatant into a pre-spun (12,000 – 16,000 x g for 1-2 min) 2 mL Phase Lock Gel tube. Add an equal volume of (Acid Phenol)/Chloroform and invert gently to mix (do not vortex).
- Spin** at 12,000-16,000 x g for 5 minutes to separate phases
- Pour** supernatant to a 1.5 mL tube. **Add** 1/10 volume 3 M NaAcetate (pH 5.2-5.6), 5 µL glycogen, and an equal volume of room temperature isopropanol. **Mix** by inversion.
- Spin** 12,000 x g at 4 °C for 20 minutes
- Rinse** pellet with 200 µL of 70% ethanol, spin 2-3 minutes. **Repeat**.
- Air dry** pellet, with tube inverted on Kimwipe, for 10 minutes
- Resuspend** pellet in ~100 µL DEPC-H₂O
- Quantitate** RNA concentrations

Make a 1:10 dilution in H₂O and measure absorbance at 260 and 280 nm (OD₂₆₀ of 1 = 40 µg/mL; 260/280 ratio should be between 1.8 and 2.0)